
Product Name: PWP1 Rabbit Monoclonal Antibody**Catalog #: AMRe84079**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 56 kDa ; Observed MW: 75 kDa

Antigen Information


Gene Name	PWP1
Alternative Names	PWP1;;PWP1
Gene ID	
SwissProt ID	Q13610
Immunogen	A synthesized peptide derived from human PWP1

Background

Chromatin-associated factor that regulates transcription.Regulates Pol I-mediated rRNA biogenesis and, probably, Pol III-mediated transcription.Regulates the epigenetic status of rDNA.

Research Area

Image Data



Western blot analysis of PWP1 expression in HepG2 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

