
Product Name: RAB45 Rabbit Monoclonal Antibody**Catalog #: AMRe84082**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100
Molecular Weight	Calculated MW: 83 kDa ; Observed MW: 85 kDa

Antigen Information


Gene Name	RAB45
Alternative Names	RAB45;;RASEF
Gene ID	
SwissProt ID	Q8IZ41
Immunogen	A synthesized peptide derived from human RASEF

Background

Binds predominantly GDP, and also GTP.Acts as a dynein adapter protein that activates dynein-mediated transport and dynein-dynactin motility on microtubules.

Research Area

Image Data



Western blot analysis of RAB45 expression in (1) A431 cell lysate; (2) Raw264.7 cell lysate.