

Product Name: PRPF4 Rabbit Monoclonal Antibody**Catalog #: AMRe84085**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	58 kDa

Antigen Information

Gene Name	PRPF4
Alternative Names	hPrp4; HPRP4P; PRP4; Prp4p; PRPF4; RP70; SNRNP60;;PRPF4
Gene ID	
SwissProt ID	O43172
Immunogen	A synthesized peptide derived from human PRPF4


Background

Plays a role in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex that is involved in spliceosome assembly,

and as component of the precatalytic spliceosome (spliceosome B complex).

Research Area

Image Data



Western blot analysis of PRPF4 expression in HeLa cell lysate.