
Product Name: PRPF4 Rabbit Monoclonal Antibody**Catalog #: AMRe84088**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	58 kDa

Antigen Information


Gene Name	PRPF4
Alternative Names	hPrp4; HPRP4P; PRP4; Prp4p; PRPF4; RP70; SNRNP60;;PRPF4
Gene ID	
SwissProt ID	O43172
Immunogen	A synthesized peptide derived from human PRPF4

Background

Plays a role in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex that is involved in spliceosome assembly, and as component of the precatalytic spliceosome (spliceosome B complex).

Research Area

Image Data



Western blot analysis of PRPF4 expression in (1) Raji cell lysate; (2) RAW 264.7 cell lysate; (3) C6 cell lysate.