
Product Name: PMF1 Rabbit Monoclonal Antibody**Catalog #: AMRe84090**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW: 23 kDa ; Observed MW: 25 kDa

Antigen Information

Gene Name	PMF1
Alternative Names	Pmf1; Polyamine modulated factor 1;;Polyamine modulated factor 1
Gene ID	
SwissProt ID	Q6P1K2
Immunogen	A synthesized peptide derived from human Polyamine modulated factor 1


Background

Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore formation during mitosis. May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of

SSAT.

Research Area

Image Data



Western blot analysis of PMF1 expression in (1) HeLa cell lysate; (2) RAW 264.7 cell lysate.