
Product Name: RPS8 Rabbit Monoclonal Antibody**Catalog #: AMRe84096**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	24 kDa

Antigen Information


Gene Name	RPS8
Alternative Names	RPS8;;RPS8
Gene ID	
SwissProt ID	P62241
Immunogen	A synthesized peptide derived from human RPS8

Background


Belongs to the ribosomal protein S8e family.

Research Area

Image Data



Western blot analysis of RPS8 expression in (1) 293 cell lysate; (2) RAW 264.7 cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)