

**Product Name: Dopamine Receptor D3 Rabbit
Monoclonal Antibody
Catalog #: AMRe84128**

Summary

Production Name	Dopamine Receptor D3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Dopamine Receptor D3
Alternative Names	D(3) dopamine receptor; Dopamine D3 receptor; DRD3; ETM1; FET1;;DRD3
Gene ID	
SwissProt ID	P35462. A synthesized peptide derived from human DRD3

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	44 kDa

Background


**Product Name: Dopamine Receptor D3 Rabbit
Monoclonal Antibody
Catalog #: AMRe84128**



Dopamine receptor whose activity is mediated by G proteins which inhibit adenylyl cyclase. Promotes cell proliferation.

Research Area

Image Data



Western blot analysis of Dopamine Receptor D3 expression in SH-SY5Y cell lysate.

Note

For research use only.