

Product Name: HOXA13 Rabbit Monoclonal Antibody**Catalog #: AMRe84131**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.61mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 40 kDa ; Observed MW: 42 kDa

Antigen Information

Gene Name	HOXA13
Alternative Names	Homeobox1J; HomeoboxA13; HOX 1; HOX 1J; HOX A13; HOX1; HOX1J; HOXA 13; HOXA13; Transcription factor HOXA13;;HOX A13
Gene ID	
SwissProt ID	P31271
Immunogen	A synthesized peptide derived from human HOX A13


Background

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific


positional identities on the anterior-posterior axis.

Research Area


Image Data



Western blot analysis of HOXA13 expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.