
Product Name: RPS19 Rabbit Monoclonal Antibody**Catalog #: AMRe84134**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
Molecular Weight	16 kDa

Antigen Information


Gene Name	RPS19
Alternative Names	DBA; DBA1; Ribosomal protein S19; rps19; S19;;RPS19
Gene ID	
SwissProt ID	P39019
Immunogen	A synthesized peptide derived from human RPS19

Background

Component of the small ribosomal subunit.The ribosome is a large ribonucleoprotein complex responsible for the synthesis of proteins in the cell.Required for pre-rRNA processing and maturation of 40S ribosomal subunits.

Research Area

Image Data



Western blot analysis of RPS19 expression in HeLa cell lysate.