

Product Name: MAT2A Rabbit Monoclonal Antibody**Catalog #: AMRe84135**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.61mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200
Molecular Weight	Calculated MW: 44 kDa ; Observed MW: 48 kDa

Antigen Information


Gene Name	MAT2A
Alternative Names	AMS2; MAT 2; MAT 2A; MAT II; Mat2a; MATA2; MATII; SAMS 2; SAMS2;;MAT2A
Gene ID	
SwissProt ID	P31153
Immunogen	A synthesized peptide derived from human MAT2A

Background


Catalyzes the formation of S-adenosylmethionine from methionine and ATP.

Research Area


Image Data



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of MAT2A expression in (1) Ramos cell lysate; (2) Mouse spleen lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.