
Product Name: Bmf Rabbit Monoclonal Antibody**Catalog #: AMRe84145**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,FC 1:20-1:100
Molecular Weight	Calculated MW: 21 kDa ; Observed MW: 24 kDa

Antigen Information


Gene Name	Bmf
Alternative Names	Bcl2 modifying factor; Bmf;;Bcl 2 modifying factor
Gene ID	
SwissProt ID	Q96LC9
Immunogen	A synthesized peptide derived from human Bcl 2 modifying factor

Background

May play a role in apoptosis. Isoform 1 seems to be the main initiator.

Research Area

Image Data



Western blot analysis of Bmf expression in Ramos cell lysate.