
Product Name: DEK Rabbit Monoclonal Antibody**Catalog #: AMRe84148**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,IP
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	43 kDa

Antigen Information


Gene Name	DEK
Alternative Names	Dek; DEK oncogene; Protein DEK;;DEK
Gene ID	
SwissProt ID	P35659
Immunogen	A synthesized peptide derived from human DEK

Background

Involved in chromatin organization.

Research Area

Image Data



Western blot analysis of DEK expression in 293T cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

