
Product Name: Nogo Receptor Rabbit Monoclonal Antibody**Catalog #: AMRe84151**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW: 51 kDa ; Observed MW: 66 kDa

Antigen Information

Gene Name	Nogo Receptor
Alternative Names	NgR; Nogo 66 receptor; Nogo receptor; Nogor; Reticulon 4 receptor; Rtn4r;;Reticulon 4 receptor
Gene ID	
SwissProt ID	Q9BZR6
Immunogen	A synthesized peptide derived from human Reticulon 4 receptor


Background

Receptor for RTN4, OMG and MAG,Functions as a receptor for the sialylated gangliosides GT1b and GM1.Besides, functions as

a receptor for chondroitin sulfate proteoglycans.

Research Area

Image Data



Western blot analysis of Nogo Receptor expression in HepG2 cell lysate.