

Product Name: CTRL Rabbit Monoclonal Antibody**Catalog #: AMRe84159**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 28 kDa ; Observed MW: 24,28 kDa

Antigen Information


Gene Name	CTRL
Alternative Names	CTRL; CTRL1;;CTRL
Gene ID	
SwissProt ID	P40313
Immunogen	A synthesized peptide derived from human CTRL

Background

This gene encodes a serine-type endopeptidase with chymotrypsin- and elastase-2-like activities. The gene encoding this zymogen is expressed specifically in the pancreas and likely functions as a digestive enzyme.

Research Area

Image Data



Western blot analysis of CTRL expression in Mouse pancreas lysate.