

---

**Product Name: CYP26A1 Rabbit Monoclonal Antibody****Catalog #: AMRe84170**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC
<b>Reactivity</b>	Human,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.61mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,IHC 1:100-1:200
<b>Molecular Weight</b>	Calculated MW: 56 kDa ; Observed MW: 50 kDa

**Antigen Information**


<b>Gene Name</b>	CYP26A1
<b>Alternative Names</b>	CP26; CYP26; cyp26a1; hP450RAI; P450RAI; P450RAI1;;Cytochrome P450 26A1
<b>Gene ID</b>	
<b>SwissProt ID</b>	O43174
<b>Immunogen</b>	A synthesized peptide derived from human Cytochrome P450 26A1

**Background**


A cytochrome P450 monooxygenase involved in the metabolism of retinoates (RAs), the active metabolites of vitamin A, and critical signaling molecules in animals.

## Research Area


## Image Data



Western blot analysis of CYP26A1 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.