
Product Name: LIMPII Rabbit Monoclonal Antibody**Catalog #: AMRe84173**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 54 kDa ; Observed MW: 77 kDa

Antigen Information


Gene Name	LIMPII
Alternative Names	AMRF; CD36; CD36L2; EPM4; HLGP85; LGP85; LIMP2; LIMPII; Scarb2; SRBII;;SCARB2
Gene ID	
SwissProt ID	Q14108
Immunogen	A synthesized peptide derived from human SCARB2

Background

Acts as a lysosomal receptor for glucosylceramidase (GBA1) targeting.

Research Area

Image Data



Western blot analysis of LIMPII expression in (1) HeLa cell lysate; (2) Mouse liver lysate.