
Product Name: NCF1 Rabbit Monoclonal Antibody**Catalog #: AMRe84203**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IP 1:20-1:50
Molecular Weight	45 kDa

Antigen Information


Gene Name	NCF1
Alternative Names	NCF1; NCF1A; NOXO2; SH3PXD1A; p47phox;;p47 phox
Gene ID	
SwissProt ID	P14598
Immunogen	A synthesized peptide derived from human p47 phox

Background

NCF2, NCF1, and a membrane bound cytochrome b558 are required for activation of the latent NADPH oxidase (necessary for superoxide production).

Research Area

Image Data



Western blot analysis of NCF1 expression in Raji cell lysate.