
Product Name: eIF2B epsilon Rabbit Monoclonal Antibody**Catalog #: AMRe84218**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
Molecular Weight	80 kDa

Antigen Information


Gene Name	eIF2B epsilon
Alternative Names	CACH; CLE; EIF 2B; EIF2B; EIF2B5; EIF2BE; EIF2Bepsilon;;eIF2B5
Gene ID	
SwissProt ID	Q13144
Immunogen	A synthesized peptide derived from human eIF2B5

Background

Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.

Research Area

Image Data



Western blot analysis of eIF2B epsilon expression in HeLa cell lysate.