

Product Name: Neuronatin Rabbit Monoclonal Antibody
Catalog #: AMRe84220



Summary

Production Name	Neuronatin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Neuronatin
Alternative Names	Neuronatin; NNAT; Peg 5; Peg5;;Neuronatin
Gene ID	
SwissProt ID	Q16517. A synthesized peptide derived from human Neuronatin

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	9 kDa

Background

May participate in the maintenance of segment identity in the hindbrain and pituitary development, and maturation or


Product Name: Neuronatin Rabbit Monoclonal Antibody
Catalog #: AMRe84220



maintenance of the overall structure of the nervous system. May function as a regulatory subunit of ion channels.

Research Area

Image Data



Western blot analysis of Neuronatin expression in SH-SY5Y cell lysate.

Note

For research use only.