

Product Name: CYP1B1 Rabbit Monoclonal Antibody**Catalog #: AMRe84252**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
Molecular Weight	61 kDa

Antigen Information


Gene Name	CYP1B1
Alternative Names	CP1B; Cyp1b1; CYP1B1; Cytochrome P450 1B1; GLC3A; P4501B1;;Cytochrome P450 1B1
Gene ID	
SwissProt ID	Q16678
Immunogen	A synthesized peptide derived from human Cytochrome P450 1B1

Background


A cytochrome P450 monooxygenase involved in the metabolism of various endogenous substrates, including fatty acids, steroid hormones and vitamins.

Research Area

Image Data



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of CYP1B1 expression in Human fetal kidney lysate.