

Product Name: DDAH2 Rabbit Monoclonal Antibody
Catalog #: AMRe84264



Summary

Production Name	DDAH2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	DDAH2
Alternative Names	DDAH; DDAH II; DDAH2; DDAHII; Dimethylargininase 2; G6a;;DDAH2
Gene ID	
SwissProt ID	O95865. A synthesized peptide derived from human DDAH2

Application

Dilution Ratio	WB:1:1000-1:2000,IP:1:20-1:50
Molecular Weight	30 kDa

Background

Putative hydrolase with unknown substrate (Probable). Does not hydrolyze N(G),N(G)-dimethyl-L-arginine (ADMA) which


Product Name: DDAH2 Rabbit Monoclonal Antibody
Catalog #: AMRe84264



acts as an inhibitor of NOS

Research Area

Image Data



Western blot analysis of DDAH2 expression in MCF7 cell lysate.

Note

For research use only.