

Summary

Production Name	HES5 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,IHC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	HES5
Alternative Names	HES 5; Hes5;;HES5
Gene ID	
SwissProt ID	Q5TA89.A synthesized peptide derived from human HES5

Application

Dilution Ratio	WB:1:1000-1:2000,IP:1:20-1:50,IHC:1:100-1:200
Molecular Weight	18 kDa

Background


Product Name: HES5 Rabbit Monoclonal Antibody
Catalog #: AMRe84266




Transcriptional repressor of genes that require a bHLH protein for their transcription. Plays an important role as neurogenesis negative regulator (By similarity).

Research Area


Image Data



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Product Name: HES5 Rabbit Monoclonal Antibody
Catalog #: AMRe84266



Note

For research use only.