

Product Name: Casein Rabbit Monoclonal Antibody**Catalog #: AMRe84267**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, IP 1:20-1:50
Molecular Weight	22 kDa

Antigen Information


Gene Name	Casein
Alternative Names	Casein;;alpha S1 casein
Gene ID	
SwissProt ID	P47710
Immunogen	A synthesized peptide derived from human alpha S1 casein

Background

Important role in the capacity of milk to transport calcium phosphate. Casoxin D acts as opioid antagonist and has vasorelaxing activity mediated by bradykinin B1 receptors.

Research Area

Image Data



Western blot analysis of Casein expression in Human milk lysate.