

**Product Name: Retinoic Acid Receptor gamma Rabbit
Monoclonal Antibody
Catalog #: AMRe84272**



Summary

Production Name	Retinoic Acid Receptor gamma Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Retinoic Acid Receptor gamma
Alternative Names	NR1B3; RAR gamma; RARC; Rarg;; RAR gamma
Gene ID	
SwissProt ID	P13631. A synthesized peptide derived from human RAR gamma

Application

Dilution Ratio	WB:1:1000-1:2000, ICC:1:50-1:200
Molecular Weight	50 kDa

Background


**Product Name: Retinoic Acid Receptor gamma Rabbit
Monoclonal Antibody
Catalog #: AMRe84272**



Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands, all-trans or 9-cis retinoic acid, and regulate gene expression in various biological processes.

Research Area

Image Data



Western blot analysis of Retinoic Acid Receptor gamma expression in HeLa cell lysate.

Note

For research use only.