

Product Name: PLTP Rabbit Monoclonal Antibody**Catalog #: AMRe84298**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.61mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 55 kDa ; Observed MW: 80 kDa

Antigen Information


Gene Name	PLTP
Alternative Names	BPIFE; HDLCQ9; OD107; Pltp;;PLTP
Gene ID	
SwissProt ID	P55058
Immunogen	A synthesized peptide derived from human PLTP

Background


Converts HDL into larger and smaller particles. May play a key role in extracellular phospholipid transport and modulation of hdl particles.

Research Area


Image Data



Western blot analysis of PLTP expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.