

Product Name: Engrailed 1 Rabbit Monoclonal Antibody
Catalog #: AMRe84306



Summary

Production Name	Engrailed 1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Engrailed 1
Alternative Names	EN1; Engrailed 1; Engrailed; Hu-En-1; Mo-en.1;;EN1
Gene ID	
SwissProt ID	Q05925.A synthesized peptide derived from human EN1

Application

Dilution Ratio	WB:1:1000-1:2000,IP:1:20-1:50,FC:1:20-1:100
Molecular Weight	Calculated MW: 40 kDa ; Observed MW: 55 kDa

Background


Product Name: Engrailed 1 Rabbit Monoclonal Antibody
Catalog #: AMRe84306




Belongs to the engrailed homeobox family.

Research Area

Image Data



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Engrailed 1 expression in 293T cell lysate.

Note

For research use only.