

Product Name: JAG2 Rabbit Monoclonal Antibody
Catalog #: AMRe84313



Summary

Production Name	JAG2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	JAG2
Alternative Names	HJ2; JAG2; Jagged2; SER2; Syndactylism;;Jagged 2
Gene ID	
SwissProt ID	Q9Y219. A synthesized peptide derived from human Jagged 2

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	Calculated MW: 133 kDa ; Observed MW: 145 kDa

Background


Putative Notch ligand involved in the mediation of Notch signaling. Involved in limb development.

Product Name: JAG2 Rabbit Monoclonal Antibody
Catalog #: AMRe84313



Research Area

Image Data



Western blot analysis of JAG2 expression in HUVEC cell lysate.

Note

For research use only.