

Summary

Production Name	BAP31 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC,FC,IF
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	BAP31
Alternative Names	BCR-associated protein Bap31; p28 Bap31; Protein CDM;;BCAP31
Gene ID	
SwissProt ID	P51572.A synthesized peptide derived from human BCAP31

Application

Dilution Ratio	WB:1:1000-1:2000,IHC:1:100-1:200;IF:1:50-1:200,ICC:1:50-1:200,FC:1:20-1:100
Molecular Weight	Calculated MW: 28 kDa ; Observed MW: 25 kDa

Background


Product Name: BAP31 Rabbit Monoclonal Antibody
Catalog #: AMRe84318




May play a role in anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi. May be involved in CASP8-mediated apoptosis.

Research Area


Image Data



Western blot analysis of BAP31 expression in NCCIT cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Product Name: BAP31 Rabbit Monoclonal Antibody
Catalog #: AMRe84318



Note

For research use only.