

Product Name: Annexin IV Rabbit Monoclonal Antibody**Catalog #: AMRe84319**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
Molecular Weight	Calculated MW: 36 kDa ; Observed MW: 30 kDa

Antigen Information


Gene Name	Annexin IV
Alternative Names	AIV; AnnexinA4; AnnexinIV; ANX 4; ANX A4; ANX4; ANXA4; Chromobindin4; EndonexinI; LipocortinIV; P32.5; P33/41; PAPII; PIG28; PP4X; Protein II; Xanx-4; ZAP36;;Annexin A4
Gene ID	
SwissProt ID	P09525
Immunogen	A synthesized peptide derived from human Annexin A4

Background

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Research Area

Image Data



Western blot analysis of Annexin IV expression in HepG2 cell lysate.