
Product Name: VEGFB Rabbit Monoclonal Antibody**Catalog #: AMRe84337**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 22 kDa ; Observed MW: 16-62 kDa

Antigen Information


Gene Name	VEGFB
Alternative Names	VEGF B; Vegfb; VEGFL; VRF;;VEGF B
Gene ID	
SwissProt ID	P49765
Immunogen	A synthesized peptide derived from human VEGF B

Background

Growth factor for endothelial cells. VEGF-B167 binds heparin and neuropilin-1 whereas the binding to neuropilin-1 of VEGF-B186 is regulated by proteolysis.

Research Area

Image Data



Western blot analysis of VEGFB expression in 293T cell transfected with VEGFB recombinant protein.