
Product Name: Pannexin 1 Rabbit Monoclonal Antibody**Catalog #: AMRe84367**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IP 1:20-1:50
Molecular Weight	48 kDa

Antigen Information

Gene Name	Pannexin 1
Alternative Names	innexin; MRS1; Pannexin 1; Panx1; PX1;;Pannexin 1
Gene ID	
SwissProt ID	Q96RD7
Immunogen	A synthesized peptide derived from human Pannexin 1


Background

The protein encoded by this gene belongs to the innexin family. Innexin family members are the structural components of gap junctions. This protein and pannexin 2 are abundantly expressed in central nerve system (CNS) and are coexpressed in various

neuronal populations.

Research Area

Image Data



Western blot analysis of Pannexin 1 expression in Caco-2 cell lysate.