

Product Name: ADAMTS13 Rabbit Monoclonal Antibody**Catalog #: AMRe84379**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, FC 1:20-1:100
Molecular Weight	154 kDa

Antigen Information


Gene Name	ADAMTS13
Alternative Names	ADAMTS13; ADAM-TS; ADAM-TS13; ADAMTS-13; ADAM-TS 13; C9orf8; TTP; VWF-CP; VWFCP;;ADAMTS13
Gene ID	
SwissProt ID	Q76LX8
Immunogen	A synthesized peptide derived from human ADAMTS13

Background

Cleaves the vWF multimers in plasma into smaller forms thereby controlling vWF-mediated platelet thrombus formation.

Research Area

Image Data



Western blot analysis of ADAMTS13 expression in A549 cell lysate.