
Product Name: MST1 Rabbit Monoclonal Antibody**Catalog #: AMRe84381**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
Molecular Weight	Calculated MW: 80 kDa ; Observed MW: 85 kDa

Antigen Information


Gene Name	MST1
Alternative Names	HGFL; MSP; MST1; NF15S2;;MST1
Gene ID	
SwissProt ID	P26927
Immunogen	A synthesized peptide derived from human MST1

Background

Probably has no proteolytic activity, since crucial characteristic of serine proteases catalytic sites are not conserved.

Research Area

Image Data



Western blot analysis of MST1 expression in HepG2 cell lysate.