

---

**Product Name: Cytochrome P450 17A1 Rabbit Monoclonal Antibody****Catalog #: AMRe84385**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ICC,FC,IP
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 57 kDa ; Observed MW: 50 kDa

**Antigen Information**


<b>Gene Name</b>	CPT7
<b>Alternative Names</b>	CPT7; CYP17; P450C17; S17AH; CYP17A1;;Cytochrome P450 17A1
<b>Gene ID</b>	
<b>SwissProt ID</b>	P05093
<b>Immunogen</b>	A synthesized peptide derived from human Cytochrome P450 17A1

**Background**

A cytochrome P450 monooxygenase involved in corticoid and androgen biosynthesis.


## Research Area

## Image Data




Western blot analysis of Cytochrome P450 17A1 expression in Jurkat cell lysate.


All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

