
Product Name: ZNF217 Rabbit Monoclonal Antibody**Catalog #: AMRe84388**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 115 kDa ; Observed MW: 125 kDa

Antigen Information


Gene Name	ZNF217
Alternative Names	ZABC1; ZN217; ZNF217;;ZNF217
Gene ID	
SwissProt ID	O75362
Immunogen	A synthesized peptide derived from human ZNF217

Background

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Research Area

Image Data



Western blot analysis of ZNF217 expression in Caco2 cell lysate.