

Product Name: MC1R eceptor Rabbit Monoclonal Antibody**Catalog #: AMRe84391**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	35 kDa

Antigen Information


Gene Name	MC1R eceptor
Alternative Names	CMM5; MC1R; MSHR; SHEP2;;MSHR
Gene ID	
SwissProt ID	Q01726
Immunogen	A synthesized peptide derived from human MSHR

Background

Receptor for MSH (alpha, beta and gamma) and ACTH. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

Research Area

Image Data



Western blot analysis of MC1 Receptor expression in A375 cell lysate.