
Product Name: Cytochrome P450 2C9 Rabbit Monoclonal Antibody**Catalog #: AMRe84424**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100
Molecular Weight	Calculated MW: 56 kDa ; Observed MW: 50 kDa

Antigen Information


Gene Name	Cytochrome P450 2C9
Alternative Names	CPC9; CYP2C; CYP2C10; CYP2C9; CYP11C9; P450 PB1; P450MP;;Cytochrome P450 2C9
Gene ID	
SwissProt ID	P11712
Immunogen	A synthesized peptide derived from human Cytochrome P450 2C9

Background

A cytochrome P450 monooxygenase involved in the metabolism of various endogenous substrates, including fatty acids and steroids.

Research Area

Image Data



Western blot analysis of Cytochrome P450 2C9 expression in HepG2 cell lysate.