

Product Name: ETV6 Rabbit Monoclonal Antibody
Catalog #: AMRe84436



Summary

Production Name	ETV6 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	ETV6
Alternative Names	ETV6; Tel; TEL oncogene; TEL/ABL; TEL1;;ETV6
Gene ID	
SwissProt ID	P41212. A synthesized peptide derived from human ETV6

Application

Dilution Ratio	WB:1:1000-1:2000, ICC:1:50-1:200, FC:1:20-1:100
Molecular Weight	53 kDa

Background

Transcriptional repressor; binds to the DNA sequence 5'-CCGGAAGT-3'. Plays a role in hematopoiesis and malignant


Product Name: ETV6 Rabbit Monoclonal Antibody
Catalog #: AMRe84436



transformation.

Research Area

Image Data



Western blot analysis of ETV6 expression in A431 cell lysate.

Note

For research use only.