

Product Name: CRABP1 Rabbit Monoclonal Antibody**Catalog #: AMRe84440**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW: 16 kDa ; Observed MW: 14 kDa

Antigen Information


Gene Name	CRABP1
Alternative Names	CRABP; Crabp1; CRABPI; RBP5;;CRABP1
Gene ID	
SwissProt ID	P29762
Immunogen	A synthesized peptide derived from human CRABP1

Background

Cytosolic CRABPs may regulate the access of retinoic acid to the nuclear retinoic acid receptors.

Research Area

Image Data



Western blot analysis of CRABP1 expression in MCF7 cell lysate.