

Product Name: PGAM1 Rabbit Monoclonal Antibody**Catalog #: AMRe84448**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,FC 1:20-1:100
Molecular Weight	29 kDa

Antigen Information

Gene Name	PGAM1
Alternative Names	PGAM1; PGAMA; PGAMB; ;PGAMA1
Gene ID	
SwissProt ID	P18669
Immunogen	A synthesized peptide derived from human PGAMA1


Background

Catalyzes the interconversion of 2-phosphoglycerate and 3-phosphoglycerate a crucial step in glycolysis, by using 2,3-bisphosphoglycerate. Also catalyzes the interconversion of (2R)-2,3-bisphosphoglycerate and (2R)-3-phospho-glyceroyl

phosphate.

Research Area

Image Data



Western blot analysis of PGAM1 expression in A431 cell lysate.