

Product Name: ESRRG Rabbit Monoclonal Antibody
Catalog #: AMRe84450



Summary

Production Name	ESRRG Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	ESRRG
Alternative Names	ERR 3; ERR gamma 2; ERR3; ERRG2; ERRgamma; ESRRG; NR3B3;;Estrogen related receptor gamma
Gene ID	
SwissProt ID	P62508.A synthesized peptide derived from human Estrogen related receptor gamma

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	51 kDa

Product Name: ESRRG Rabbit Monoclonal Antibody
Catalog #: AMRe84450




Background


Orphan receptor that acts as a transcription activator in the absence of bound ligand. Binds specifically to an estrogen response element and activates reporter genes controlled by estrogen response elements.

Research Area

Image Data



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Estrogen Related Receptor gamma expression in Human fetal heart lysate.

Note

For research use only.