
Product Name: Placental lactogen Rabbit Monoclonal Antibody**Catalog #: AMRe84454**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	25 kDa

Antigen Information

Gene Name	Placental lactogen
Alternative Names	CS-2; CS1; CSA; CSB; CSH1; CSH2; CSMT; hCS-B; hCSA; Lactogen; PL; PL1; Placental lactogen;;Placental lactogen
Gene ID	
SwissProt ID	P01243
Immunogen	A synthesized peptide derived from human Placental lactogen


Background

Produced only during pregnancy and is involved in stimulating lactation, fetal growth and metabolism. Does not interact with

GHR but only activates PRLR through zinc-induced dimerization.

Research Area

Image Data



Western blot analysis of Placental lactogen expression in JAR cell lysate.