

Product Name: COL11A1 Rabbit Monoclonal Antibody
Catalog #: AMRe84457



Summary

Production Name	COL11A1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	COL11A1
Alternative Names	COL11A1; COLL6; Collagen alpha 1; collagen XI alpha 1; STL2; STL3;;Collagen XI alpha 1
Gene ID	
SwissProt ID	P12107. A synthesized peptide derived from human Collagen XI alpha 1

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	Calculated MW: 181 kDa ; Observed MW: 150 kDa

Background


May play an important role in fibrillogenesis by controlling lateral growth of collagen II fibrils.

Product Name: COL11A1 Rabbit Monoclonal Antibody
Catalog #: AMRe84457



Research Area

Image Data



Western blot analysis of COL11A1 expression in K562 cell lysate.

Note

For research use only.