

Product Name: KHSRP Rabbit Monoclonal Antibody**Catalog #: AMRe84500**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 73 kDa ; Observed MW: 82 kDa

Antigen Information

Gene Name	KHSRP
Alternative Names	FBP2; FUBP2; KHSRP; KSRP; p75;;KHSRP
Gene ID	
SwissProt ID	Q92945
Immunogen	A synthesized peptide derived from human KHSRP


Background

Binds to the dendritic targeting element and may play a role in mRNA trafficking.Part of a ternary complex that binds to the downstream control sequence (DCS) of the pre-mRNA. Mediates exon inclusion in transcripts that are subject to tissue-specific

alternative splicing. May interact with single-stranded DNA from the far-upstream element (FUSE).

Research Area

Image Data



Western blot analysis of KHSRP expression in HeLa cell lysate.