
Product Name: PROS1 Rabbit Monoclonal Antibody**Catalog #: AMRe84504**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	75 kDa

Antigen Information

Gene Name	PROS1
Alternative Names	Preproprotein S; PROS; proS1; PS21; PS22; PS23; PS24; PS25; PS26; PSA; THPH5; THPH6;;PROS1
Gene ID	
SwissProt ID	P07225
Immunogen	A synthesized peptide derived from human PROS1


Background

Anticoagulant plasma protein; it is a cofactor to activated protein C in the degradation of coagulation factors Va and VIIIa. It

helps to prevent coagulation and stimulating fibrinolysis.

Research Area

Image Data



Western blot analysis of PROS1 expression in Human plasma lysate.