

Product Name: TPO Rabbit Monoclonal Antibody**Catalog #: AMRe84508**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
Molecular Weight	103 kDa

Antigen Information


Gene Name	TPO
Alternative Names	MSA; TDH2A; TPO; TPX;;Thyroid peroxidase
Gene ID	
SwissProt ID	P07202
Immunogen	A synthesized peptide derived from human Thyroid peroxidase

Background

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T3 and T4.

Research Area

Image Data



Western blot analysis of TPO expression in Human thyroid lysate.