

Summary

Production Name	SCGN Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,IHC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	SCGN
Alternative Names	CALBL; Scgn; SECRET; Secretagogin; SEGN; Setagin;;Secretagogin
Gene ID	
SwissProt ID	O76038.A synthesized peptide derived from human Secretagogin

Application

Dilution Ratio	WB:1:1000-1:2000,IP:1:20-1:50,IHC:1:100-1:200
Molecular Weight	32 kDa

Background


Product Name: SCGN Rabbit Monoclonal Antibody
Catalog #: AMRe84509




The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation.

Research Area


Image Data



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Product Name: SCGN Rabbit Monoclonal Antibody
Catalog #: AMRe84509



Note

For research use only.