

Product Name: BCL2L15 Rabbit Monoclonal Antibody**Catalog #: AMRe84520**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.61mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
Molecular Weight	Calculated MW: 18 kDa ; Observed MW: 52 kDa

Antigen Information

Gene Name	BCL2L15
Alternative Names	Bcl-2-like protein 15; Bcl2-L-15; BCL2-like 15; Bcl2l15; Bflk; Gm566; Pro apoptotic Bcl 2 protein; ;Bcl 2 L15
Gene ID	
SwissProt ID	Q5TBC7
Immunogen	A synthesized peptide derived from human Bcl 2 L15


Background

Predicted to be involved in apoptotic process and regulation of apoptotic process. Predicted to be active in cytosol and


nucleus.

Research Area


Image Data



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of BCL2L15 in His-tagged BCL2L15 cell lysate.