
Product Name: FYB Rabbit Monoclonal Antibody**Catalog #: AMRe84523**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
Molecular Weight	Calculated MW: 85 kDa ; Observed MW: 120 kDa

Antigen Information


Gene Name	FYB
Alternative Names	ADAP; Fyb; FYN binding protein; p120/p130; SLAP130;;FYN binding protein 1
Gene ID	
SwissProt ID	O15117
Immunogen	A synthesized peptide derived from human FYN binding protein 1

Background

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T-cells.May play a role in linking T-cell signaling to remodeling of the actin cytoskeleton.

Research Area

Image Data



Western blot analysis of FYB expression in TF-1 cell lysate.