
Product Name: GSTK1 Rabbit Monoclonal Antibody**Catalog #: AMRe84544**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	25 kDa

Antigen Information

Gene Name	GSTK1
Alternative Names	GST; GST class kappa; GST13; Gstk1; hGSTK1;;GSTK1
Gene ID	
SwissProt ID	Q9Y2Q3
Immunogen	A synthesized peptide derived from human GSTK1


Background

Glutathione S-transferase that catalyzes the conjugation of glutathione to exogenous and endogenous compounds. Significant

glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene.

Research Area

Image Data



Western blot analysis of GSTK1 expression in 293T cell lysate.